



Serving Scientific Community  
Since 1985



**TECHNO PHARMCHEM**

An ISO 9001 : 2008 Certified Company

Registered Off. : 101-A, Deep Enclave, Pocket "D", Ashok Vihar, Phase-III, Delhi-110 052 (India)  
Corporate Off. : 152, Vardhman City Centre, Near Shakti Nagar Underbridge, Delhi-110052 (India)  
Works : Plot No.1022, Modern Industrial Estate, Bahadurgarh-124507, Haryana (India)  
Telephone : 0091 - 11 - 23646422 E-mail : tecpharm@gmail.com  
Website : www.technopharmchem.com

## Nutrient Broth

---

### Product Code : CM0149

Nutrient Broth is used for the cultivation of less fastidious microorganisms, can be enriched with blood or other biological fluids.

### Composition\*\*

Ingredients	Gms / Litre
Peptone Bacteriological	5.000
Sodium chloride	5.000
Meat Extract	1.500
Yeast extract	1.500

Final pH ( at 25°C) 7.4±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 13 grams in 1000 ml distilled water. Heat, if necessary, to dissolve the medium completely. Dispense as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

### Principle And Interpretation

Nutrient media are basic culture media used for maintaining microorganisms, cultivating fastidious organisms by enriching with serum or blood and are also used for purity checking prior to biochemical or serological testing. Nutrient Agar is ideal for demonstration and teaching purposes where a more prolonged survival of cultures at ambient temperature is often required without risk of overgrowth that can occur with more nutritious substrate. This relatively simple formula has been retained and is still widely used in the microbiological examination of variety of materials and is also recommended by standard methods. It is one of the several non-selective media useful in routine cultivation of microorganisms. It can be used for the cultivation and enumeration of bacteria which are not particularly fastidious. Addition of different biological fluids such as horse or sheep blood, serum, egg yolk etc. makes it suitable for the cultivation of related fastidious organisms. Peptone, Meat extract and yeast extract provide the necessary nitrogen compounds, carbon, vitamins and also some trace ingredients necessary for the growth of bacteria. Sodium chloride maintains the osmotic equilibrium of the medium.

### Quality Control

Appearance	:	Cream to yellow homogeneous free flowing powder
Colour and Clarity of prepared medium	:	Light amber colored clear solution in tubes
Reaction	:	Reaction of 1.3% w/v aqueous solution at 25°C. pH : 7.4±0.2
pH	:	7.20-7.60



Serving Scientific Community  
Since 1985



**TECHNO PHARMCHEM**

**An ISO 9001 : 2008 Certified Company**

Registered Off. : 101-A, Deep Enclave, Pocket "D", Ashok Vihar, Phase-III, Delhi-110 052 (India)  
Corporate Off. : 152, Vardhman City Centre, Near Shakti Nagar Underbridge, Delhi-110052 (India)  
Works : Plot No.1022, Modern Industrial Estate, Bahadurgarh-124507, Haryana (India)  
Telephone : 0091 - 11 - 23646422 E-mail : tecpharm@gmail.com  
Website : www.technopharmchem.com

### Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

Cultural Response

Organism	Inoculum (CFU)	Growth
Escherichia coli ATCC 25922	50-100	good-luxuriant
Pseudomonas aeruginosa ATCC 27853	50-100	good-luxuriant
Staphylococcus aureus ATCC 25923	50-100	good-luxuriant
Streptococcus pyogenes ATCC 19615	50-100	good-luxuriant

### Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C.  
Use before expiry date on the label.